

VG500

The INFICON Vacuum Gauge Display VGD500 connects to our CDG product line to provide a 4 Digit Display. It's a small unit which Displays the pressure of the gauge. The pressure reading is selectable in pressure unit Torr, mTorr, mbar and Pa.



利点

- 4 Digit display for easy read out
- Large 10mm active LED display - readable from distance and wide range gauge of angle
- In-line plug design
- 15pin D-Sub, no extra power connection
- Compact small size
- Compatible with all Fullscale ≥ 100 mTorr
- pressure Units selectable (Torr, mTorr, mbar and Pa)
- CE certified

注文の情報

タイプ

VG500 Vacuum Gauge Display

VG500

399-653

仕様

タイプ		VGD500
Display		4 digits
Connection		
Gauge Side		D-Sub, 15-pin, female
Measurement cable side		D-Sub, 15-pin, male
Signal		digital input RS232
Pressure unit (adjustable)		Torr (default), mTorr, mbar, Pa
Supply		
Power consumption	W	≤1
Supply		
Voltage		+14 ...+30V (dc) or +/-15V (ac)
Degree of protection		IP 40
Admissible temperature		
Operation		+5 ... +50°C
Storage		-20 ... +85°C
Use		indoors only
Connectable gauges in the measurement range (F.S.) 0.1mbar/Torr to 1000mbar/Torr		CDG025D / -S, CDG025D-X3, CDG045D ... CDG200D, CDG045D2 ... CDG100D2
Weight		58g
Dimensions		50x63x34mm



www.inficon.com reachus@inficon.com

製品向上のプログラムを続けるため、仕様は通知なく変更されます。
RateWatchはINFICONの商標です。その他すべての商標は、それぞれの所有者の資産となります。

(2016-08) © 2016 INFICON